

For Research Use Only

KRT19 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86250-8-PBS



Basic Information

Catalog Number:

86250-8-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

BC002539

GeneID (NCBI):

3880

UNIPROT ID:

P08727

Full Name:

keratin 19

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

250894D5

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02308-2, Human Cytokeratin 19 Recombinant Matched Antibody Pair - PBS only. 86250-8-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard SY02555. 86250-7-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL



The mean Cytokeratin 19 concentration was determined to be 228.05 ng/mL in MCF-7 cell extract based on a 1.30 mg/mL extract load.