

For Research Use Only

# LBP Recombinant antibody

Catalog Number: 86253-1-RR



## Basic Information

Catalog Number: 86253-1-RR	GenBank Accession Number: NM_008489.2	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 16803	CloneNo.: 250583C3
Source: Rabbit	UNIPROT ID: Q61805	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: lipopolysaccharide binding protein	
Immunogen Catalog Number: EG3163	Calculated MW: 53 kDa	
	Observed MW: 60 kDa	

## Applications

Tested Applications: WB, ELISA	Positive Controls:
Species Specificity: mouse	WB : mouse serum, mouse ovary tissue

## Background Information

LBP (lipopolysaccharide-binding protein) is an acute-phase protein, predominantly synthesized by hepatocytes. LBP plays an essential role in the body's innate immune response, whereby it protects against infections by recognizing and binding bacterial lipopolysaccharides to activate immune cells.

## Storage

**Storage:**  
Store at -20°C. Stable for one year after shipment.  
**Storage Buffer:**  
PBS with 0.02% sodium azide and 50% glycerol, pH7.3  
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)

W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86253-1-RR (Lbp antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.