

For Research Use Only

# RPS3A Recombinant antibody, PBS Only (Detector)

Catalog Number: 86254-1-PBS



## Basic Information

Catalog Number:

86254-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG5273

GenBank Accession Number:

BC066926

GeneID (NCBI):

6189

UNIPROT ID:

P61247

Full Name:

ribosomal protein S3A

Calculated MW:

30 kDa

Purification Method:

Protein A purification

CloneNo.:

251033G5

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean RPS3A concentration was determined to be 48.1 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02339-1, RPS3A Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86254-2-PBS. Detection antibody: 86254-1-PBS. Standard: Ag5273. Range: 0.625-40 ng/mL.