

For Research Use Only

CNPY3 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86256-1-PBS



Basic Information

Catalog Number:

86256-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7382

GenBank Accession Number:

BC004423

GeneID (NCBI):

10695

UNIPROT ID:

Q9BT09

Full Name:

canopy 3 homolog (zebrafish)

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

250998F12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CNPY3 concentration was determined to be 20.6 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02342-1, CNPY3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86256-2-PBS. Detection antibody: 86256-1-PBS. Standard: Ag7382. Range: 0.469-60 ng/mL.