

For Research Use Only

# JAKMIP1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86257-1-PBS



## Basic Information

Catalog Number:

86257-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4816

GenBank Accession Number:

BC047075

GeneID (NCBI):

152789

UNIPROT ID:

Q96N16

Full Name:

janus kinase and microtubule interacting protein 1

Calculated MW:

73 kDa

Purification Method:

Protein A purification

CloneNo.:

251026D1

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean JAKMIP1 concentration was determined to be 76.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02332-1, JAKMIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86257-2-PBS. Detection antibody: 86257-1-PBS. Standard: Ag4816. Range: 0.313-20 ng/mL.