

For Research Use Only

JAKMIP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86257-2-PBS



Basic Information

Catalog Number:

86257-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG4816

GenBank Accession Number:

BC047075

GeneID (NCBI):

152789

UNIPROT ID:

Q96N16

Full Name:

janus kinase and microtubule
interacting protein 1

Calculated MW:

73 kDa

Purification Method:

Protein A purification

CloneNo.:

251026F1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean JAKMIP1 concentration was determined to be 76.9 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02332-1, JAKMIP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86257-2-PBS. Detection antibody: 86257-1-PBS. Standard: Ag4816. Range: 0.313-20 ng/mL.