

For Research Use Only

C22orf28 Recombinant antibody, PBS Only (Capture)

Catalog Number:86258-3-PBS



Basic Information

Catalog Number: 86258-3-PBS	GenBank Accession Number: BC000151	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 51493	CloneNo.: 250866B12
Source: Rabbit	UNIPROT ID: Q9Y3I0	
Isotype: IgG	Full Name: chromosome 22 open reading frame 28	
Immunogen Catalog Number: AG13836	Calculated MW: 55 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02340-2, Human RTCB Recombinant Matched Antibody Pair - PBS only. 86258-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag13836. 86258-1-PBS was HRP conjugated as the detection antibody. Range: 0.391-25 ng/mL



The mean RTCB concentration was determined to be 78.30 ng/mL in HeLa cell extract based on a 1.1 mg/mL extract load and 77.70 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.