

For Research Use Only

SNRPA Recombinant antibody, PBS Only (Detector)

Catalog Number: 86259-1-PBS



Basic Information

Catalog Number:

86259-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0283

GenBank Accession Number:

BC000405

GeneID (NCBI):

6626

UNIPROT ID:

P09012

Full Name:

small nuclear ribonucleoprotein

polypeptide A

Calculated MW:

31 kDa

Purification Method:

Protein A purification

CloneNo.:

251055G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean SNRPA concentration was determined to be 81.7 ng/mL in Hela cell extract based on a 0.5 mg/mL extract load.



Cytometric bead array standard curve of MP02344-1, SNRPA Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86259-2-PBS. Detection antibody: 86259-1-PBS. Standard: Ag0283. Range: 0.781-100 ng/mL