

For Research Use Only

UMPS Recombinant antibody, PBS Only (Detector)

Catalog Number: 86260-1-PBS



Basic Information

Catalog Number:

86260-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG6620

GenBank Accession Number:

BC000364

GeneID (NCBI):

7372

UNIPROT ID:

P11172

Full Name:

uridine monophosphate synthetase

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

251020B6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean UMPS concentration was determined to be 34.7 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02341-1, UMPS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86260-2-PBS. Detection antibody: 86260-1-PBS. Standard: Ag6620. Range: 0.469-60 ng/mL.