

For Research Use Only

CXCL3/GRO Gamma Recombinant antibody, PBS Only

Catalog Number: 86261-1-PBS



Basic Information

Catalog Number:

86261-1-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2982

GenBank Accession Number:

NM_002090.3

GeneID (NCBI):

2921

UNIPROT ID:

P19876

Full Name:

chemokine (C-X-C motif) ligand 3

Calculated MW:

11 kDa

Observed MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

250114C9

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human, mouse, rat

Background Information

CXC chemokine ligand 3 (CXCL3) is a member of CXC-type chemokine family that is identified as a major regulator in immune and inflammation responses. CXCL3 is broadly expressed in various human tumor types, and it is also known to play a critical role in mediating tumor development and progression.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 86261-1-RR (CXCL3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86261-1-PBS in a different storage buffer formulation.