

For Research Use Only

CXCL3/GRO Gamma Recombinant antibody

Catalog Number: 86261-1-RR



Basic Information

Catalog Number: 86261-1-RR	GenBank Accession Number: NM_002090.3	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 2921	CloneNo.: 250114C9
Source: Rabbit	UNIPROT ID: P19876	Recommended Dilutions: WB: 1:5000-1:50000
Isotype: IgG	Full Name: chemokine (C-X-C motif) ligand 3	
Immunogen Catalog Number: EG2982	Calculated MW: 11 kDa	
	Observed MW: 11 kDa	

Applications

Tested Applications: WB, ELISA	Positive Controls: WB : RAW 264.7 cells, A549 cells, rat lung tissue
Species Specificity: human, mouse, rat	

Background Information

CXC chemokine ligand 3 (CXCL3) is a member of CXC-type chemokine family that is identified as a major regulator in immune and inflammation responses. CXCL3 is broadly expressed in various human tumor types, and it is also known to play a critical role in mediating tumor development and progression.

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



RAW 264.7 cells were subjected to SDS PAGE followed by western blot with 86261-1-RR (CXCL3 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.