

For Research Use Only

DR4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86267-1-PBS



Basic Information

Catalog Number:

86267-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4600

GenBank Accession Number:

BC012866

GeneID (NCBI):

8797

UNIPROT ID:

O00220

Full Name:

tumor necrosis factor receptor
superfamily, member 10a

Calculated MW:

50 kDa

Purification Method:

Protein A purification

CloneNo.:

250737D3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DR4 concentration was determined to be 9.4 ng/mL in Hela cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02323-1, DR4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86267-2-PBS. Detection antibody: 86267-1-PBS. Standard: Eg4600. Range: 0.781-100 ng/mL.