

For Research Use Only

Lp-PLA2/PLA2G7 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86271-1-PBS



Basic Information

Catalog Number:

86271-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3756

GenBank Accession Number:

NM_005084.4

GeneID (NCBI):

7941

UNIPROT ID:

Q13093

Full Name:

phospholipase A2, group VII (platelet-activating factor acetylhydrolase, plasma)

Calculated MW:

50kDa

Purification Method:

Protein A purification

CloneNo.:

250253D12

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Lp-PLA2 concentration was determined to be 1615.8 ng/mL in human serum.



Cytometric bead array standard curve of MP02319-1, Lp-PLA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86271-2-PBS. Detection antibody: 86271-1-PBS. Standard: Eg3756. Range: 0.781-100 ng/mL.