

For Research Use Only

FNDC3B Recombinant antibody, PBS Only (Capture)

Catalog Number: 86272-2-PBS



Basic Information

Catalog Number:

86272-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG18357

GenBank Accession Number:

BC039297

GeneID (NCBI):

64778

UNIPROT ID:

Q53EP0

Full Name:

fibronectin type III domain containing 3B

Calculated MW:

1204 aa, 133 kDa

Purification Method:

Protein A purification

CloneNo.:

250827F2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FNDC3B concentration was determined to be 70.0 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02305-1, FNDC3B Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86272-2-PBS. Detection antibody: 86272-1-PBS. Standard: Ag18357. Range: 3.125-200 ng/mL