

For Research Use Only

# Transthyretin Recombinant antibody, PBS Only (Capture)

Catalog Number: 86273-2-PBS



## Basic Information

Catalog Number:

86273-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG0931

GenBank Accession Number:

NM\_000371.4

GeneID (NCBI):

7276

UNIPROT ID:

P02766

Full Name:

transthyretin

Calculated MW:

16 kDa

Purification Method:

Protein A purification

CloneNo.:

250820E7

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Prealbumin/TTR concentration was determined to be 172.8 ug/mL in human serum.



Cytometric bead array standard curve of MP02304-1, Prealbumin/TTR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86273-2-PBS. Detection antibody: 86273-1-PBS. Standard: Eg0931. Range: 1.56-50 ng/mL