

For Research Use Only

Neurocan Recombinant antibody, PBS Only

Catalog Number: 86279-1-PBS



Basic Information

Catalog Number:

86279-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4217

GenBank Accession Number:

NM_004386.2

GeneID (NCBI):

1463

UNIPROT ID:

O14594

Full Name:

neurocan

Calculated MW:

143 kDa

Observed MW:

270 kDa

Purification Method:

Protein A purification

CloneNo.:

250695D11

Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



SH-SY5Y cells were subjected to SDS PAGE followed by western blot with 86279-1-RR (Neurocan antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86279-1-PBS in a different storage buffer formulation.



fetal human brain tissue were subjected to SDS PAGE followed by western blot with 86279-1-RR (Neurocan antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86279-1-PBS in a different storage buffer formulation.