

For Research Use Only

# Complement C6 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86286-1-PBS



## Basic Information

Catalog Number:

86286-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG4529

GenBank Accession Number:

BC035723

GeneID (NCBI):

729

UNIPROT ID:

P13671

Full Name:

complement component 6

Calculated MW:

934 aa, 105 kDa

Purification Method:

Protein A purification

CloneNo.:

250971G12

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Complement C6 concentration was determined to be 17.3 ug/mL in human serum.



Cytometric bead array standard curve of MP02310-1, Complement C6 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86286-2-PBS. Detection antibody: 86286-1-PBS. Standard: Eg4529. Range: 0.313-40 ng/mL.