

For Research Use Only

HDGF2 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86289-2-PBS



Basic Information

Catalog Number:

86289-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7281

GenBank Accession Number:

BC000755

GeneID (NCBI):

84717

UNIPROT ID:

Q7Z4V5

Full Name:

hepatoma-derived growth factor-related protein 2

Calculated MW:

74 kDa

Purification Method:

Protein A purification

CloneNo.:

250810D6

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HDGF2 concentration was determined to be 13.5 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02306-1, HDGF2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86289-2-PBS. Detection antibody: 86289-1-PBS. Standard: Ag7281. Range: 0.781-100 ng/mL.