

For Research Use Only

# PNN Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86292-3-PBS



## Basic Information

Catalog Number:

86292-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG38940

GenBank Accession Number:

BC062602

GeneID (NCBI):

5411

UNIPROT ID:

Q9H307

Full Name:

pinin, desmosome associated protein

Calculated MW:

82 kDa

Purification Method:

Protein A purification

CloneNo.:

251028E4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



The mean Pinin concentration was determined to be 8.3 ng/mL in HepG2 cell extract based on a 0.9 mg/mL extract load.



Cytometric bead array standard curve of MP02380-1, Pinin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86292-3-PBS. Detection antibody: 86292-2-PBS. Standard: Ag38940. Range: 0.156-10 ng/mL.