

For Research Use Only

TOX Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86295-1-PBS



Basic Information

Catalog Number:

86295-1-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3171

GenBank Accession Number:

BC016665

GeneID (NCBI):

9760

UNIPROT ID:

O94900

Full Name:

thymocyte selection-associated high
mobility group box

Calculated MW:

526 aa, 58 kDa

Purification Method:

Protein A purification

CloneNo.:

251078B3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean TOX concentration was determined to be 2016.4 ng/mL in EL-4 cell extract based on a 3.7 mg/mL extract load.



Cytometric bead array standard curve of MP02373-1, TOX Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86295-2-PBS. Detection antibody: 86295-1-PBS. Standard: Ag3171. Range: 0.781-100 ng/mL