

For Research Use Only

IL13RA2,CD213A2 Recombinant antibody, PBS Only (Detector)

Catalog Number:86298-1-PBS



Basic Information

Catalog Number:

86298-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2579

GenBank Accession Number:

BC020739

GeneID (NCBI):

3598

UNIPROT ID:

Q14627

Full Name:

interleukin 13 receptor, alpha 2

Calculated MW:

380 aa, 44 kDa

Purification Method:

Protein A purification

CloneNo.:

250952G5

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IL-13RA2 concentration was determined to be 75.2 ng/mL in A375 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02369-1, IL-13RA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86298-2-PBS. Detection antibody: 86298-1-PBS. Standard: Eg2579. Range: 0.781-100 ng/mL.