

For Research Use Only

LCP1 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86299-2-PBS



Basic Information

Catalog Number: 86299-2-PBS	GenBank Accession Number: BC010271	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 3936	CloneNo.: 251063E1
Source: Rabbit	UNIPROT ID: P13796	
Isotype: IgG	Full Name: lymphocyte cytosolic protein 1 (L-plastin)	
Immunogen Catalog Number: AG3571	Calculated MW: 627 aa, 64 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean LCP1 concentration was determined to be 4.5 ng/mL in K-562 cell extract based on a 3.0 mg/mL extract load.



Cytometric bead array standard curve of MP02368-1, LCP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86299-2-PBS. Detection antibody: 86299-1-PBS. Standard: Ag3571. Range: 0.469-60 ng/mL.