

For Research Use Only

# PRPF31 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86315-3-PBS



## Basic Information

Catalog Number:

86315-3-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG26716

GenBank Accession Number:

NM\_015629

GeneID (NCBI):

26121

UNIPROT ID:

Q8WWY3

Full Name:

PRP31 pre-mRNA processing factor 31 homolog (S. cerevisiae)

Calculated MW:

55 kDa

Purification Method:

Protein A purification

CloneNo.:

251005B10

## Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean PRPF31 concentration was determined to be 14.26 ng/mL in HeLa cell extract based on a 1.10 mg/mL extract load.



Sandwich ELISA standard curve of MP02586-1, Human PRPF31 Recombinant Matched Antibody Pair - PBS only. 86315-1-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag26716. 86315-3-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL