

For Research Use Only

Mouse CD39 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86327-3-PBS



Basic Information

Catalog Number:

86327-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3483

GenBank Accession Number:

NM_009848.4

GeneID (NCBI):

12495

UNIPROT ID:

P55772

Full Name:

ectonucleoside triphosphate
diphosphohydrolase 1

Calculated MW:

57kDa

Purification Method:

Protein A purification

CloneNo.:

250764C7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD39/ENTPD1 concentration was determined to be 263.3 ng/mL in Mouse spleen tissue extract based on a 6.4 mg/mL extract load.



Cytometric bead array standard curve of MP02378-1, MOUSE CD39/ENTPD1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86327-3-PBS. Detection antibody: 86327-2-PBS. Standard: Eg3483. Range: 0.469-60 ng/mL.