

For Research Use Only

Mouse CD39 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86327-4-PBS



Basic Information

Catalog Number:

86327-4-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3483

GenBank Accession Number:

NM_009848.4

GeneID (NCBI):

12495

UNIPROT ID:

P55772

Full Name:

ectonucleoside triphosphate
diphosphohydrolase 1

Calculated MW:

57 kDa

Purification Method:

Protein A purification

CloneNo.:

25076486

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02378-2, Mouse CD39/ENTPD1 Recombinant Matched Antibody Pair - PBS only. 86327-4-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg3483. 86327-1-PBS was HRP conjugated as the detection antibody. Range: 78.1-5000 pg/mL.



The mean CD39/ENTPD1 concentration was determined to be 926.3 ng/mL in mouse spleen extract based on a 8.5 mg/mL extract load and 2,212.6 ng/mL in mouse lung extract based on a 5.9 mg/mL extract load.