

For Research Use Only

QARS Recombinant antibody, PBS Only (Detector)

Catalog Number: 86333-1-PBS



Basic Information

Catalog Number:

86333-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3375

GenBank Accession Number:

BC016634

GeneID (NCBI):

5859

UNIPROT ID:

P47897

Full Name:

glutaminyI-tRNA synthetase

Calculated MW:

775 aa, 88 kDa

Purification Method:

Protein A purification

CloneNo.:

251052C4

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GlnRS concentration was determined to be 55.2 ng/mL in PC-3 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02363-1, GlnRS Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86333-2-PBS. Detection antibody: 86333-1-PBS. Standard: Ag3375. Range: 0.781-100 ng/mL.