

For Research Use Only

QARS Recombinant monoclonal antibody

Catalog Number:86333-3-RR



Basic Information

Catalog Number: 86333-3-RR	GenBank Accession Number: BC016634	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 5859	CloneNo.: 251052A2
Isotype: IgG	UNIPROT ID: P47897	Recommended Dilutions: WB: 1:5000-1:50000 IF/ICC: 1:200-1:800
Immunogen Catalog Number: AG3375	Full Name: glutaminyl-tRNA synthetase	
	Calculated MW: 775 aa, 88 kDa	
	Observed MW: 88 kDa	

Applications

Tested Applications: WB, IF/ICC, ELISA	Positive Controls:
Species Specificity: human, mouse, rat	WB : HeLa cells, HepG2 cells, PC-3 cells, Ramos cells, Raji cells, NIH/3T3 cells, HSC-T6 cells
	IF/ICC : NIH/3T3 cells,

Background Information

GlnRS (glutaminyl-tRNA synthetase) is a class I aminoacyl-tRNA synthetase. Aminoacyl-tRNA synthetases are enzymes that charge tRNAs with their cognate amino acids. The specificity of this reaction determines the fidelity of mRNA translation. At least 1 synthetase exists in the cytoplasm for each amino acid. GlnRS is present in eukaryotes but is absent from many prokaryotes, mitochondria, and chloroplasts, in which Gln-tRNA(Gln) is formed by transamidation of the misacylated Glu-tRNA(Gln)(PMID:8078941).

Storage

Storage:
Store at -20°C. Stable for one year after shipment.
Storage Buffer:
PBS with 0.02% sodium azide and 50% glycerol, pH7.3
Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86333-3-RR (QARS antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Immunofluorescent analysis of (-20°C Methanol) fixed NIH/3T3 cells using QARS antibody (86333-3-RR, Clone: 251052A2) at dilution of 1:400 and CoraLite® 488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2).



Biolayer interferometry (BLI) kinetic assays of 86333-3-RR against Human QARS were performed. The affinity constant is 0.176 nM.