

For Research Use Only

# HSD17B4 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86334-1-PBS



## Basic Information

Catalog Number:

86334-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7165

GenBank Accession Number:

BC003098

GeneID (NCBI):

3295

UNIPROT ID:

P51659

Full Name:

hydroxysteroid (17-beta)  
dehydrogenase 4

Calculated MW:

80 kDa

Purification Method:

Protein A purification

CloneNo.:

251058C4

## Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

## Background Information

### Storage

Storage:

Store at -80°C.

**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**

## Selected Validation Data



The mean HSD17B4 concentration was determined to be 3562.2 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02361-1, HSD17B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86334-2-PBS. Detection antibody: 86334-1-PBS. Standard: Ag7165. Range: 1.562-200 ng/mL.