

For Research Use Only

HSD17B4 Recombinant antibody, PBS Only (Capture)

Catalog Number: 86334-2-PBS



Basic Information

Catalog Number:

86334-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7165

GenBank Accession Number:

BC003098

GeneID (NCBI):

3295

UNIPROT ID:

P51659

Full Name:

hydroxysteroid (17-beta)
dehydrogenase 4

Calculated MW:

80 kDa

Purification Method:

Protein A purification

CloneNo.:

251058E1

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean HSD17B4 concentration was determined to be 3562.2 ng/mL in HeLa cell extract based on a 1.5 mg/mL extract load.



Cytometric bead array standard curve of MP02361-1, HSD17B4 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86334-2-PBS. Detection antibody: 86334-1-PBS. Standard: Ag7165. Range: 1.562-200 ng/mL.