

For Research Use Only

Mouse LAMP-1/CD107a Recombinant antibody, PBS Only (Detector)

Catalog Number: 86337-1-PBS



Basic Information

Catalog Number:

86337-1-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3858

GenBank Accession Number:

NM_010684.3

GeneID (NCBI):

16783

UNIPROT ID:

P11438

Full Name:

lysosomal-associated membrane protein 1

Calculated MW:

44 kDa

Purification Method:

Protein A purification

CloneNo.:

250972H2

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CD107a / LAMP1 concentration was determined to be 28.9 ng/mL in NIH/3T3 cell extract based on a 1.4 mg/mL extract load.



Cytometric bead array standard curve of MP02362-1, MOUSE CD107a / LAMP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86337-2-PBS. Detection antibody: 86337-1-PBS. Standard: Eg3858. Range: 0.625-80 ng/mL.