

For Research Use Only

α -Actinin Recombinant antibody, PBS Only (Detector)

Catalog Number: 86340-1-PBS



Basic Information

Catalog Number:	86340-1-PBS	GenBank Accession Number:	BC015766	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	87	CloneNo.:	250833B7
Source:	Rabbit	UNIPROT ID:	P12814		
Isotype:	IgG	Full Name:	actinin, alpha 1		
Immunogen Catalog Number:	AG1859	Calculated MW:	103 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Alpha Actinin concentration was determined to be 254.6 ng/mL in NIH/3T3 cell extract based on a 1.4 mg/mL extract load.



Cytometric bead array standard curve of MP02360-1, Alpha Actinin Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86340-2-PBS. Detection antibody: 86340-1-PBS. Standard: Ag1859. Range: 0.781-100 ng/mL.