

For Research Use Only

TMEM66 Recombinant antibody, PBS Only (Capture)

Catalog Number:86350-3-PBS



Basic Information

Catalog Number: 86350-3-PBS	GenBank Accession Number: BC015012	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 51669	CloneNo.: 250805C3
Source: Rabbit	UNIPROT ID: Q96BY9	
Isotype: IgG	Full Name: transmembrane protein 66	
Immunogen Catalog Number: AG29388	Calculated MW: 339 aa, 37 kDa	

Applications

Tested Applications:
Sandwich ELISA, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Sandwich ELISA standard curve of MP02452-1, Human TMEM66 Recombinant Matched Antibody Pair - PBS only. 86350-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag29388. 86350-2-PBS was HRP conjugated as the detection antibody. Range: 1.56-100 ng/mL



The mean TMEM66 concentration was determined to be 314.83 ng/mL in A549 cell extract based on a 1.1 mg/mL