

For Research Use Only

PJA2 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86366-4-PBS



Basic Information

Catalog Number:

86366-4-PBS

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG3434

GenBank Accession Number:

BC030826

GeneID (NCBI):

9867

UNIPROT ID:

O43164

Full Name:

praja ring finger 2

Calculated MW:

708 aa, 78 kDa

Purification Method:

Protein A purification

CloneNo.:

251077A2

Applications

Tested Applications:

Sandwich ELISA, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean PJA2 concentration was determined to be 2.40 ng/mL in SW480 cell extract based on a 3.0 mg/mL extract load and 4.47 ng/mL in Jurkat cell extract based on a 1.2 mg/mL extract load.



Sandwich ELISA standard curve of MP02550-1, Human PJA2 Recombinant Matched Antibody Pair - PBS only. 86366-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Ag3434. 86366-4-PBS was HRP conjugated as the detection antibody. Range: 0.195-12.5 ng/mL