

For Research Use Only

# Mouse PR3/PRTN3 Recombinant antibody, PBS Only

Catalog Number: 86384-3-PBS



## Basic Information

Catalog Number:

86384-3-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2750

GenBank Accession Number:

NM\_011178.2

GeneID (NCBI):

19152

UNIPROT ID:

Q61096

Full Name:

proteinase 3

Calculated MW:

28 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

250371B12

## Applications

Tested Applications:

WB, Indirect ELISA

Species Specificity:

mouse, rat

## Background Information

### Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

## Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 86384-3-RR (Prtn3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86384-3-PBS in a different storage buffer formulation.



rat spleen tissue were subjected to SDS PAGE followed by western blot with 86384-3-RR (Prtn3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86384-3-PBS in a different storage buffer formulation.