

For Research Use Only

PR3/PRTN3 Recombinant antibody

Catalog Number: 86384-3-RR



Basic Information

Catalog Number:

86384-3-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2750

GenBank Accession Number:

NM_011178.2

GeneID (NCBI):

19152

UNIPROT ID:

Q61096

Full Name:

proteinase 3

Calculated MW:

28 kDa

Observed MW:

28 kDa

Purification Method:

Protein A purification

CloneNo.:

250371B12

Recommended Dilutions:

WB: 1:500-1:2000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

mouse, rat

Positive Controls:

WB : mouse spleen tissue, rat spleen tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



mouse spleen tissue were subjected to SDS PAGE followed by western blot with 86384-3-RR (Prtn3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.



rat spleen tissue were subjected to SDS PAGE followed by western blot with 86384-3-RR (Prtn3 antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours.