

For Research Use Only

Mouse Igfbp1 Recombinant antibody, PBS Only (Capture)

Catalog Number:86385-2-PBS



Basic Information

Catalog Number: 86385-2-PBS	GenBank Accession Number: NM_008341.4	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 16006	CloneNo.: 250310B2
Source: Rabbit	UNIPROT ID: P47876	
Isotype: IgG	Full Name: insulin-like growth factor binding protein 1	
Immunogen Catalog Number: EG2934	Calculated MW: 30 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
mouse

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean IGFBP1 concentration was determined to be 66.2 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02401-1, MOUSE IGFBP1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86385-2-PBS. Detection antibody: 86385-1-PBS. Standard: Eg2934. Range: 0.312-40 ng/mL.