

For Research Use Only

IL-25/IL-17E Recombinant monoclonal antibody

Catalog Number: 86396-1-RR



Basic Information

Catalog Number:

86396-1-RR

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG3316

GenBank Accession Number:

NM_080729.3

GeneID (NCBI):

140806

UNIPROT ID:

Q8VHH8

Full Name:

interleukin 25

Calculated MW:

19KD

Observed MW:

20 kda

Purification Method:

Protein A purification

CloneNo.:

250033G5

Recommended Dilutions:

WB: 1:1000-1:4000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human, mouse, rat

Positive Controls:

WB : HT-29 cells, rat testis tissue, LoVo cells, SW480 cells, mouse colon tissue, mouse testis tissue, rat colon tissue

Background Information

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86396-1-RR (IL-25/IL-17E antibody) at dilution of 1:2000 incubated at room temperature for 1.5 hours.



rat testis tissue were subjected to SDS PAGE followed by western blot with 86396-1-RR (IL-25/IL-17E antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours.