

For Research Use Only

Phospho-Tyrosine Recombinant monoclonal antibody, PBS Only

Catalog Number: 86398-1-PBS



Basic Information

Catalog Number:

86398-1-PBS

Source:

Rabbit

Isotype:

IgG

GenBank Accession Number:

GeneID (NCBI):

Full Name:

Purification Method:

Protein A purification

CloneNo.:

250495G7

Applications

Tested Applications:

WB, IF/ICC, Indirect ELISA

Species Specificity:

-

Background Information

Phospho-Tyrosine antibody refers to a class of pan-specific antibodies or detection reagents that recognize all phosphorylated tyrosine residues regardless of the surrounding amino-acid sequence. These tools are widely used to monitor global tyrosine phosphorylation events in signal-transduction studies, cancer research and drug-discovery screens.

This antibody does not cross-react with peptides containing Phospho-Serine or Phospho-Threonine.

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.


Selected Validation Data



Non-treated A431 cells and EGF (HZ-1326) treated A431 cells were subjected to SDS PAGE followed by western blot with 86398-1-RR (Phospho-Tyrosine antibody) at dilution of 1:1000 incubated at room temperature for 1.5 hours. The membrane was stripped and re-blotted with GAPDH (60004-1-Ig) antibody as a loading control. This data was developed using the same antibody clone with 86398-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed hEGF treated A431 cells using Phospho-Tyrosine antibody (86398-1-RR, Clone: 250495G7) at dilution of 1:100 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86398-1-PBS in a different storage buffer formulation.



Immunofluorescent analysis of (4% PFA) fixed hEGF treated A431 cells using Phospho-Tyrosine antibody (86398-1-RR, Clone: 250495G7) at dilution of 1:200 and CoraLite®488-Conjugated Goat Anti-Rabbit IgG(H+L) (SA00013-2). This data was developed using the same antibody clone with 86398-1-PBS in a different storage buffer formulation.