

For Research Use Only

MIA3 Recombinant monoclonal antibody, PBS Only (Capture)

Catalog Number:86403-2-PBS



Basic Information

Catalog Number: 86403-2-PBS	GenBank Accession Number: BC047116	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 375056	CloneNo.: 250698B9
Isotype: IgG	UNIPROT ID: Q5JRA6	
Immunogen Catalog Number: AG39740	Full Name: melanoma inhibitory activity family, member 3	
	Calculated MW: 214 kDa	

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean MIA3 concentration was determined to be 82.5 ng/mL in HeLa cell extract based on a 1.2 mg/mL extract load.



Cytometric bead array standard curve of MP02435-1, MIA3 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86403-2-PBS. Detection antibody: 86403-1-PBS. Standard: Ag39740. Range: 0.625-20 ng/mL.