

For Research Use Only

Mouse TYRO3 Recombinant antibody, PBS Only

Catalog Number: 86404-3-PBS



Basic Information

Catalog Number: 86404-3-PBS	GenBank Accession Number: NM_019392.2	Purification Method: Protein A purification
Concentration: 1 mg/ml	GeneID (NCBI): 22174	CloneNo.: 250330812
Source: Rabbit	UNIPROT ID: P55144-1	
Isotype: IgG	Full Name: TYRO3 protein tyrosine kinase 3	
Immunogen Catalog Number: EG2918	Calculated MW: 96 kDa	
	Observed MW: 130 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
mouse, rat

Background Information

TYRO3 has a structure includes a cytoplasmic domain and an extracellular ligand binding region. And ligand binding at the cell surface induces dimerization and autophosphorylation of TYRO3 on its intracellular domain that provides docking sites for downstream signaling molecules. TYRO3 is associated with many physiological processes including cell survival, migration and differentiation. It has been reported that TYRO3 can be expressed in various regions of the mammalian central nervous system. The molecular mass of TYRO3 is 130 kDa, the apparent molecular mass is higher due to glycosylation. (PMID: 8873131, 23354874, 21539887)

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C
Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86404-3-RR (Tyro3 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86404-3-PBS in a different storage buffer formulation.