

For Research Use Only

# PLXNC1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number:86406-2-PBS



## Basic Information

<b>Catalog Number:</b> 86406-2-PBS	<b>GenBank Accession Number:</b> NM_005761.3	<b>Purification Method:</b> Protein A purification
<b>Source:</b> Rabbit	<b>GeneID (NCBI):</b> 10154	<b>CloneNo.:</b> 250374D2
<b>Isotype:</b> IgG	<b>UNIPROT ID:</b> O60486	
<b>Immunogen Catalog Number:</b> EG2721	<b>Full Name:</b> plexin C1	
	<b>Calculated MW:</b> 176 kDa	

## Applications

**Tested Applications:**  
Sandwich ELISA, Indirect ELISA

**Species Specificity:**  
human

## Background Information

### Storage

**Storage:**  
Store at -80°C.  
**The product is shipped with ice packs. Upon receipt, store it immediately at -80°C**

**Storage Buffer:**  
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: [Proteintech-CN@ptglab.com](mailto:Proteintech-CN@ptglab.com)


W: [ptgcn.com](http://ptgcn.com)

**This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.**


## Selected Validation Data



The mean PLXNC1 concentration was determined to be 2.30 ng/mL in HepG2 cell extract based on a 1.20 mg/mL extract load.



Sandwich ELISA standard curve of MP02634-1, Human PLXNC1 Recombinant Matched Antibody Pair - PBS only. 86406-3-PBS was coated to a plate as the capture antibody and incubated with serial dilutions of standard Eg2721. 86406-2-PBS was HRP conjugated as the detection antibody. Range: 0.156-10 ng/mL.



Bi-layer interferometry (BLI) kinetic assay of 86406-2-PBS against Human PLXNC1 was performed. The affinity constant is 1.38 nM.