

For Research Use Only

DLK1 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86410-2-PBS



Basic Information

Catalog Number:	86410-2-PBS	GenBank Accession Number:	NM_003836.6	Purification Method:	Protein A purification
Concentration:	1 mg/ml	GeneID (NCBI):	8788	CloneNo.:	250840A7
Source:	Rabbit	UNIPROT ID:	P80370-1		
Isotype:	IgG	Full Name:	delta-like 1 homolog (Drosophila)		
Immunogen Catalog Number:	EG2909	Calculated MW:	41 kDa		

Applications

Tested Applications:
Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean DLK1 concentration was determined to be 45.1 ng/mL in human serum.



Cytometric bead array standard curve of MP02426-1, DLK1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86410-3-PBS. Detection antibody: 86410-2-PBS. Standard: Eg2909. Range: 0.625-40 ng/mL.