

For Research Use Only

Casein Kinase 1 Delta Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86423-1-PBS



Basic Information

Catalog Number: 86423-1-PBS	GenBank Accession Number: BC003558	Purification Method: Protein A purification
Source: Rabbit	GeneID (NCBI): 1453	CloneNo.: 251238G5
Isotype: IgG	UNIPROT ID: P48730	
Immunogen Catalog Number: AG5856	Full Name: casein kinase 1, delta	
	Calculated MW: 47 kDa	

Applications

Tested Applications:
IF/ICC, Cytometric bead array, Indirect ELISA

Species Specificity:
human

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean CSNK1D concentration was determined to be 12.5 ng/mL in HEK-293 cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP02443-1, CSNK1D Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86423-3-PBS. Detection antibody: 86423-1-PBS. Standard: Ag5856. Range: 0.625-80 ng/mL.