

For Research Use Only

SLFN5 Recombinant antibody

Catalog Number: 86424-2-RR



Basic Information

Catalog Number:

86424-2-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG9959

GenBank Accession Number:

BC021238

GeneID (NCBI):

162394

UNIPROT ID:

Q08AF3

Full Name:

schlafen family member 5

Calculated MW:

338 aa, 38 kDa

Observed MW:

90-100 kDa

Purification Method:

Protein A purification

CloneNo.:

251214B9

Recommended Dilutions:

WB: 1:5000-1:50000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : T-47D cells, SKOV-3 cells, HaCaT cells, MDA-MB-231 cells, LNCaP cells

Background Information

Schlafen family member 5 (SLFN5) is a member of the Schlafen family of proteins, a group of type 1 IFN-inducible proteins. In addition to an AAA (ATPase) domain and a specific SLFN box, the SLFN5 gene contains a helicase domain as well as a nuclear localization sequence, which suggests a role for this protein in transcription-related processes. A growing number of findings suggest that human Schlafen family member 5 (SLFN5) plays a key role in malignancy.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86424-2-RR (SLFN5 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.



Biolayer interferometry (BLI) kinetic assays of 86424-2-RR against Human SLFN5 were performed. The affinity constant is below 1 pM.