

For Research Use Only

FGL1 Recombinant monoclonal antibody, PBS Only (Detector)

Catalog Number: 86426-3-PBS



Basic Information

Catalog Number:

86426-3-PBS

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2765

GenBank Accession Number:

NM_004467.4

GeneID (NCBI):

2267

UNIPROT ID:

Q08830

Full Name:

fibrinogen-like 1

Calculated MW:

36 kDa

Purification Method:

Protein A purification

CloneNo.:

251192G7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FGL1 concentration was determined to be 69.0 ng/mL in human plasma.



Cytometric bead array standard curve of MP02439-1, FGL1 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86426-4-PBS. Detection antibody: 86426-3-PBS. Standard: Eg2765. Range: 0.625-80 ng/mL.