

For Research Use Only

GSTA2 Recombinant antibody, PBS Only (Detector)

Catalog Number: 86428-2-PBS



Basic Information

Catalog Number:

86428-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG0781

GenBank Accession Number:

BC002895

GeneID (NCBI):

2939

UNIPROT ID:

P09210

Full Name:

glutathione S-transferase alpha 2

Calculated MW:

26 kDa

Purification Method:

Protein A purification

CloneNo.:

251223E7

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean GSTA2 concentration was determined to be 135.2 ng/mL in HepG2 cell extract based on a 1.1 mg/mL extract load.



Cytometric bead array standard curve of MP02440-1, GSTA2 Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86428-3-PBS. Detection antibody: 86428-2-PBS. Standard: Ag0781. Range: 0.781-50 ng/mL.