

For Research Use Only

FDXR Recombinant antibody, PBS Only (Capture)

Catalog Number: 86431-3-PBS



Basic Information

Catalog Number:

86431-3-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

AG7145

GenBank Accession Number:

BC063493

GeneID (NCBI):

2232

UNIPROT ID:

P22570

Full Name:

ferredoxin reductase

Calculated MW:

54 kDa

Purification Method:

Protein A purification

CloneNo.:

251211G11

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

human

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean FDXR concentration was determined to be 14.2 ng/mL in U-251 cell extract based on a 1.6 mg/mL extract load.



Cytometric bead array standard curve of MP02464-1, FDXR Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86431-3-PBS. Detection antibody: 86431-2-PBS. Standard: Ag7145. Range: 1.562-200 ng/mL.