

For Research Use Only

UFSP2 Recombinant antibody, PBS Only

Catalog Number: 86436-1-PBS



Basic Information

Catalog Number: 86436-1-PBS	GenBank Accession Number: BC010493	Purification Method: Protein A purification
Concentration: 1000 µg/ml	GeneID (NCBI): 55325	CloneNo.: 251182F9
Source: Rabbit	UNIPROT ID: Q9NUQ7	
Isotype: IgG	Full Name: UFM1-specific peptidase 2	
Immunogen Catalog Number: AG10692	Calculated MW: 469 aa, 53 kDa	
	Observed MW: 53 kDa	

Applications

Tested Applications:
WB, Indirect ELISA

Species Specificity:
human, mouse, pig

Background Information

Storage

Storage:
Store at -80°C.
The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:
PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



Various lysates were subjected to SDS PAGE followed by western blot with 86436-1-RR (UFSP2 antibody) at dilution of 1:5000 incubated at room temperature for 1.5 hours. This data was developed using the same antibody clone with 86436-1-PBS in a different storage buffer formulation.



Biolayer interferometry (BLI) kinetic assays of 86436-1-RR against Human UFSP2 were performed. The affinity constant is below 1 pM.