

For Research Use Only

XCL1 Recombinant antibody

Catalog Number: 86437-1-RR



Basic Information

Catalog Number:

86437-1-RR

Concentration:

1000 µg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2933

GenBank Accession Number:

NM_002995.2

GeneID (NCBI):

6375

UNIPROT ID:

P47992

Full Name:

chemokine (C motif) ligand 1

Calculated MW:

13 kDa

Observed MW:

11 kDa

Purification Method:

Protein A purification

CloneNo.:

250232D1

Recommended Dilutions:

WB: 1:5000-1:10000

Applications

Tested Applications:

WB, ELISA

Species Specificity:

human

Positive Controls:

WB : PMA, Ionomycin and Brefeldin A treated NK-92 cells,

Background Information

XCL1 also known as Lymphotactin, is a member of the chemokine superfamily. Chemokines function in inflammatory and immunological responses, inducing leukocyte migration and activation.

Storage

Storage:

Store at -20°C. Stable for one year after shipment.

Storage Buffer:

PBS with 0.02% sodium azide and 50% glycerol, pH7.3

Aliquoting is unnecessary for -20°C storage

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com

W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



PMA, Ionomycin and Brefeldin A treated NK-92 cells were subjected to SDS PAGE followed by western blot with 86437-1-RR (XCL1 antibody) at dilution of 1:10000 incubated at room temperature for 1.5 hours.