

For Research Use Only

Mouse Chga Recombinant antibody, PBS Only (Detector)

Catalog Number: 86447-2-PBS



Basic Information

Catalog Number:

86447-2-PBS

Concentration:

1 mg/ml

Source:

Rabbit

Isotype:

IgG

Immunogen Catalog Number:

EG2760

GenBank Accession Number:

NM_007693.2

GeneID (NCBI):

12652

UNIPROT ID:

P26339

Full Name:

chromogranin A

Calculated MW:

52 kDa

Purification Method:

Protein A purification

CloneNo.:

250318F3

Applications

Tested Applications:

Cytometric bead array, Indirect ELISA, Sample test

Species Specificity:

mouse

Background Information

Storage

Storage:

Store at -80°C.

The product is shipped with ice packs. Upon receipt, store it immediately at -80°C

Storage Buffer:

PBS only, pH7.3

For technical support and original validation data for this product please contact:


T: 4006900926

E: Proteintech-CN@ptglab.com


W: ptgcn.com

This product is exclusively available under Proteintech Group brand and is not available to purchase from any other manufacturer.

Selected Validation Data



The mean Chromogranin a concentration was determined to be 77.6 ng/mL in mouse serum.



Cytometric bead array standard curve of MP02444-1, MOUSE Chromogranin a Recombinant Matched Antibody Pair, PBS Only. Capture antibody: 86447-1-PBS. Detection antibody: 86447-2-PBS. Standard: Eg2760. Range: 0.625-80 ng/mL.